

CHEMICAL LISTINGS

These are some of the chemicals which are retained by our GP general purpose filter. When ordering a cabinet, please indicate which chemicals/organics/inorganics you are using, so that we can determine the best filter combination for your purpose.

Acetaldehyde	Cyclohexanol	Hexylene (1-hexene)	Octalene
Acetic acid	Cyclohexanone	Hexyne (butyl acetylene)	Octane
Acetic anhydride	Cyclohexene	Hydrogen selenide	Organic chemicals
Acetone			Organic vapours
Acetonitrile	Decane	Indole	Ozone
Acids (organic)	Dibromoethane	Inorganic vapours	
Acrolein	Dichlorobenzene	Iodine	Pentane
Acrylic acid	Dichlorodifluoromethane	Iodoform	Pentanol
Acrylonitrile	Dichloroethane	Iso-butane	Pentanone
Adhesives	Dichloroethylene	Isocyanates	Pentene
Alcohol	Dichloromonofluoromethane	Isoprene	Perchloroethylene
Aldehydes	Dichloronitroethane	Isopropyl acetate	Phenol
Aliphatic hydrocarbons	Dichloropropane	Isopropyl alcohol	Phosgene
Allyl alcohol	Dichlorotetrafluoroethane		Plasticisers
Amyl acetate	Diesel fumes	Kerosene	Propane
Amyl alcohol	Diethyl amine	Keytones	Propanol
Amylene	Diethyl ketone		Propionaldehyde
Aniline	Dimethylaniline	Lactic acid	Propyl acetate
Aromatic hydrocarbons	Dimethylsulphate		Propyl alcohol
Ashphalt fumes	Dioexane	Mesityl oxide	Propyl chloride
Automobile exhaust fumes	Dipentane	Methanol	Propylene
	Dipropylketone	Methyl acetate	Pyridine
Benzene		Methyl acrylate	
Bromine	Epichlorhydrin	Methylal	Rubber odour
Butadiene	Essential oils	Methyl alcohol	
Butane	Esters	Methyl bromide	Skatole
Butanol	Ethanol	Methyl butyl ketone	Solvents
Butyl acetate	Ethyl acetate	Methyl cellosolve	Styrene monomer
Butyl alcohol	Ethyl acrylate	Methyl chloride	Sulphur trioxide
Butyl cellosolve	Ethyl alcohol	Methyl chloroform	
Butyl chloride	Ethyl amine	Methyl cyclohexane	Tar fumes
Butyraldehyde	Ethyl benzene	Methyl cyclohexanol	Tetrachloroethane
Butyric acid	Ethyl bromide	Methylcyclohexanone	Tetrachloroethylene
	Ethyl chloride	Methylene chloride	Tetrahydrofuran
Camphor	Ethylene chlorohydrin	Methyl ethyl ketone	Toluene
Caprylic acid	Ethylene dichloride	Methyl formate	Toluene di isocyanate
Carbon disulphide	Ethylene oxide	Methyl isobutyl ketone	Toluidine
Carbon tetrachloride	Ethyl formate	Methyl isocyanate	Trichlorobenzene
Carbonyl sulphide	Ethyl silicate	Methyl methacrylate	Trichloroethane
Cellosolve (ethyl glycol)	Eucalyptole	Monochlorobenzene	Trichloroethylene
Chlorine			Trichlorofluoromethane
Chlorobenzene	Fluorotrichloromethane	Naphthalene	Triethylamine
Chlorobutadiene	Formaldehyde	Nicotine	Trimethylbenzene
Chloroethene	Freons	Nitro benzene	Turpentine
Chloroform		Nitroethane	
Chloronitropropane	Gasoline	Nitrogen compounds	Ureauric acid
Chloropicrin	Glutaraldehyde	Nitroglycerine	
Chloropropane		Nitromethane	Valeric acid
Corrosive gases	Halogens	Nitropropane	Valeraldehyde
Creosote	Heptane	Nitrotoluene	Vinyl acetate
Cresol	Heptanol	Nonane	Vinyl chloride
Crontonaldehyde	Heptene		Volatile organic compounds
Cumene	Hexanol		
Cyclohexane	Hexene		Xylene



Felcon Limited

Felcon Limited & Goldsworth Medical Unit 11 - 11B Euro Business Park New Road Newhaven East Sussex BN9 0DQ England
Telephone: +44 (0)1273 513434 Email: labsales@felcon.co.uk www.felcon.co.uk

Felcon has a continuous development programme aimed at further improving and developing its products. All specifications are, therefore, subject to change.
Unauthorised use is prohibited. E&OE. All trade marks acknowledged. © 2017 Felcon Limited.