



The Felcon Educational Cabinet 'Schools 800' has been designed as a mobile, easy to clean cabinet, with good all-round visibility enabling the user to clearly demonstrate experiments within the classroom environment.

- ◆ Robust simple design for use within the classroom
- ◆ Mobile, fitted with foot master castors
- ◆ Swan neck tap / sink / gas tap and auxillary electrical socket
- ◆ F/GPS carbon filtration for the removal of commonly used fumes and vapours used for general experiments in the classroom environment
- ◆ All clear acrylic panels for good all round visibility
- ◆ Easy to clean, corrosion free, all polypropylene construction
- ◆ Simple to operate

These cabinets can also be supplied as bench mounted.

Pre-filter to 5 microns in size; 95% efficient.

Multi-purpose F/GPS carbon filter for the effective removal of most commonly used solvent and acid vapours.



## 'SCHOOLS 800' EDUCATIONAL CABINET

Dimensions (W x L x H): External: 800 x 610 x 1690mm. Internal: 710 x 548 x 610mm  
Face Velocity 0.5m/s



Felcon Limited

Felcon Limited & Goldsworth Medical Unit 11 - 11B Euro Business Park New Road Newhaven East Sussex BN9 0DQ England  
Telephone: +44 (0)1273 513434 Email: labsales@felcon.co.uk www.felcon.co.uk

Felcon has a continuous development programme aimed at further improving and developing its products. All specifications are, therefore, subject to change.  
Unauthorised use is prohibited. E&OE. All trade marks acknowledged. © 2017 Felcon Limited.