

## POWDER WEIGHING STATIONS

The Felcon Powder Weighing Stations come in two sizes; 900mm wide and 1200mm wide, designed with dynamic airflow patterns, providing the best operator protection whilst also maintaining balance accuracy.

### Features:

- ◆ Designed with operator protection in mind
- ◆ Robust construction, floor mounted, ergonomic design
- ◆ Dynamic aerofoils (stainless steel)
- ◆ Solid black base and rear baffle design
- ◆ Easy to operate
- ◆ LED lighting
- ◆ Quiet fan less than 55db
- ◆ Easy to maintain and easy to clean
- ◆ Face velocity 0.3m/s – 0.5m/s
- ◆ Pre-filter to 5 microns in size
- ◆ HEPA filter to 0.5 microns in size – 99.995% efficient
- ◆ Low airflow alarm
- ◆ Airflow monitor
- ◆ Good access and all round visibility
- ◆ Auxillary electrical socket
- ◆ Safe waste disposal
- ◆ Robust stainless steel stand with pull out shelf/ work surface and foot master castors

All of these cabinets can be ducted to the outside, or can re-circulate into the laboratory.



### POWDER WEIGHING STATIONS

Dimensions (W x L x H):

External dimensions: 900mm x 780mm x 1654mm high  
1200mm x 780mm x 1654mm high

Internal dimensions: 880mm x 580mm x 730mm  
1180mm x 580mm x 730mm



Felcon Limited

Felcon Limited & Goldsworth Medical Unit 11 - 11B Euro Business Park New Road Newhaven East Sussex BN9 0DQ England  
Telephone: +44 (0)1273 513434 Email: labsales@felcon.co.uk www.felcon.co.uk

Felcon has a continuous development programme aimed at further improving and developing its products. All specifications are, therefore, subject to change.  
Unauthorised use is prohibited. E&OE. All trade marks acknowledged. © 2017 Felcon Limited.