

Bringing colour into your laboratories

Using proven carbon filter technology, Felcon introduces you to their new range of 'Sirocco' re-circulating fume hoods, offered in three standard sizes; 900mm, 1200mm and 1500mm wide. These are all polypropylene and acrylic construction, so easy to clean with no corrosion. Of course, if you don't want the colour, these can be manufactured in all clear acrylic, and offer a whole host of features as standard.

- ◆ Operator protection, both carbon and HEPA filtration
- ◆ Ergonomic design, robust construction
- ◆ Trolley mounted, with under bench storage and foot master castors
- ◆ LED lighting
- ◆ Quiet fan; less than 55db
- ◆ Dynamic air flow patterns, dynamic stainless steel aerofoil
- ◆ Airflow monitor
- ◆ Radiused edges, easy to clean
- ◆ Easy to maintain
- ◆ Non-corrosive
- ◆ 0.5m/s airflow

Custom sizes available on request.

FELCON Sirocco



FELLAB RE-CIRCULATING FUME HOODS



FELLAB RE-CIRCULATING FUME HOODS

Dimensions (W x L x H):

Sirocco 900	External dimensions: 900mm x 610mm x 1660mm Internal dimensions: 830mm x 535mm x 745mm
Sirocco 1200	External dimensions: 1200mm x 610mm x 1660mm Internal dimensions: 1130mm x 535mm x 745mm
Sirocco 1500	External dimensions: 1500mm x 610mm x 1660mm Internal dimensions: 1430mm x 535mm x 745mm



Felcon Limited

Felcon Limited & Goldsworth Medical Unit 11 - 11B Euro Business Park New Road Newhaven East Sussex BN9 0DQ England
Telephone: +44 (0)1273 513434 Email: labsales@felcon.co.uk www.felcon.co.uk

Felcon has a continuous development programme aimed at further improving and developing its products. All specifications are, therefore, subject to change.
Unauthorised use is prohibited. E&OE. All trade marks acknowledged. © 2017 Felcon Limited.