

ROBOTICS ENCLOSURES

Felcon have been enclosing automated processes within the pharmaceutical manufacturing industry since the 1970's, so it should come as no surprise that we are experts in this field and now make enclosures for laboratory robotics, with successful installations throughout Europe and the USA.

Knowing that each installation is unique, we will listen to your requirements, and carefully project manage every stage of your needs with close attention to detail.

Robotics Enclosures designed to enclose your automation, providing:

- ◆ *Product protection*
- ◆ *Operator protection*
- ◆ *Clean air environment*
- ◆ *Environmental control*
- ◆ *Custom built to meet your exact requirements*
- ◆ *Easy all round access for routine accessibility and maintenance*
- ◆ *LED lighting*
- ◆ *Colour touchscreen control panel*
- ◆ *Low fan noise*
- ◆ *Robust construction*
- ◆ *Easily maintained and serviced, readily available filters*
- ◆ *Enclosures can either be filtered or ducted to your exhaust system*



Accessories

Electrical Lift Trolley (docking station)

- ◆ Battery operated variable height trolley
- ◆ Maximum 50–300 Kgs lift
- ◆ 50 lifts per charge
- ◆ Silent operation
- ◆ Small footprint for ease of manoeuvrability
- ◆ Digital display with 3 memory settings available
- ◆ Spare wall mountable rechargeable battery packs available enabling a fully charged battery to always be ready for use
- ◆ Battery condition indicator and low alarm



Under bench polypropylene mobile storage cabinets

A range of under bench polypropylene mobile storage cabinets, sized to suit your applications. Easy to clean, and providing excellent access to equipment and reagents.



Height adjustable stainless steel trollies/ benches

Felcon manufacture a complete range of height adjustable, robust, stainless steel benches, ideal for holding and moving heavy equipment, such as liquid handlers and incubators. Sizes available on request.



Felcon Limited

Felcon Limited & Goldsworth Medical Unit 11 - 11B Euro Business Park New Road Newhaven East Sussex BN9 0DQ England
Telephone: +44 (0)1273 513434 Email: labsales@felcon.co.uk www.felcon.co.uk

Felcon has a continuous development programme aimed at further improving and developing its products. All specifications are, therefore, subject to change.
Unauthorised use is prohibited. E&OE. All trade marks acknowledged. © 2017 Felcon Limited.